

SECTION	2	RECEIPT & INSPECTION
	2.1	The Basic Vixen Kit
	2.2	Optional Items
	2.3	Inspection & Corrective Action

2 RECEIPT & INSPECTION

2.1 The Basic Vixen Kit

The Vixen wheelchair is delivered in packaging intended to protect the product during transport. Batteries where supplied will be packed separately.

The basic Vixen kit includes the following items:

1	x	Seat Unit
2	x	Arm rests, one fitted with the controller mounting, where applicable
1	x	Chassis Unit
1	x	Battery Compartment cover
2	x	Battery connector leads
1	x	Pair of footrest assemblies
1	x	User Guide

2.2 Optional Items

Depending on the model and options, one or more of the following items may also be included:

1	x	Controller
2	x	Batteries (packed separately)
1	x	Battery Charger

Headrest, Tray and any other accessories, will be supplied as ordered.

2.3 Inspection & Corrective Action

Open the cartons, remove the foam packaging and check the contents against the Despatch Advice, and check for any damage that may have occurred in transit.

Important

Any shortages must be notified to Newton Products Ltd immediately.

Transit Damage:

In accordance with our Conditions of Sale the carrier must be informed of transit damage within 48 hours of receipt of goods.